


Product Name: Recombinant Human UBE2Z Protein
Catalog #: PEH4611



Summary

Name	UBE2Z_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H832
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WDGERTAPQC LLRIKRDIMS IYKEPPPGMF VVPDVTDMTK IHALITGPFD TPYEGGFFLF VFRCPPDYPI HPPRVKLMTT GNNTVRFNPN FYRNGKVCLS ILGTWTGPAW SPAQSISSVL ISIQSLMTEN PYHNEPGFEQ ERHPGDSKNY NECIRHETIR VAVCDMMEGK CPCPEPLRGV MEKSFLEYD FYEVACKDRL HLQGQTMQDP FGKRGHFDY QSLLMRLGLI RQKVLERLHN ENAEMDSDSS

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 Z (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme Z) (Uba6-specific E2 conjugating enzyme 1) (Use1) (Ubiquitin carrier protein Z) (Ubiquitin-protein ligase Z)
--------------------------	---

Product Name: Recombinant Human UBE2Z Protein
Catalog #: PEH4611



Note

For Research Use Only , Not for Diagnostic Use