


Product Name: Recombinant Human FGF20 Protein
Catalog #: PEH4613



Summary

Name	FGF20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9NP95
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGGLEGLGQQ VGSFLLPPA GERPLLGER RSAAERSARG GPGAAQLAHL HGILRRRQLY CRTGFHLQIL PDGSVQGTRQ DHSFGILEF ISVAVGLVSI RGVDSGLYLG MNDKGELYGS EKLTSCEFIR EQFEENWYNT YSSNIYKHGD TGRRYFVALN KDGTPRDGAR SKRHQKFTHF LRPVDPERV PELYKDLLMY

image



Background

Alternative Names Fibroblast growth factor 20 (FGF-20)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGF20 Protein
Catalog #: PEH4613

