

Product Name: Recombinant Human ACP6 Protein
Catalog #: PEH4615



Summary

Name	PPA6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9NPH0
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGQLTKVGMQ QMFALGERLR KNYVEDIPFL SPTFNPQEVF IRSTNIFRNL ESTRCLLAGL FQCQKEGP II IHTDEADSEV LYPNYQSCWS LRQRTRGRRQ TASLQPGISE DLKKVKDRMG IDSSDKVDFF ILLDNVAAEQ AHNLPSCPML KRFARMIEQR AVDTSLYILP KEDRESLQMA VGPFLHILES NLLKAMDSAT APDKIRKLYL YAAHDVTFIP LLMTLGIFDH KWPPFAVDLT MELYQHLESK

image



Background

Alternative Names	Lysophosphatidic acid phosphatase type 6 (EC 3.1.3.2) (Acid phosphatase 6, lysophosphatidic) (Acid phosphatase-like protein 1) (PACPL1)
--------------------------	---

Note

Product Name: Recombinant Human ACP6 Protein
Catalog #: PEH4615



For Research Use Only , Not for Diagnostic Use