


Product Name: Recombinant Human FANCF Protein
Catalog #: PEH4617



Summary

Name	FANCF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9NPI8
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRHIHRRFGR HGPIRTALER RLHNQWRQEG GFGRGPVPL ANFQALGHCD VLLSLRLEN RALGDAARYH LVQQLFPGPG VRDADEETLQ ESLARLARRR SAVHMLRFNG YRENPNLQED SLMKTQAELL LERLQEVGKA EAERPARFLS SLWERLPQNN FLKVIALL QPPLSRRPQE ELEPGIHKSP GEGSQVLVHW LLGNSEVFAA FCRALPAGLL TLVTSRHPAL SPVYLGLLTD WGQRLHYDLQ KGIWVGTESQ DVPWEELHNR

image



Background

Alternative Names Fanconi anemia group F protein (Protein FAF)

Product Name: Recombinant Human FANCF Protein
Catalog #: PEH4617



Note

For Research Use Only , Not for Diagnostic Use