

Summary

Name	CRAC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQ79
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTQLNYGVAV TDVDHDGDFE IVVAGYNGPN LVLKYDRAQK RLVNIAVDER SSPYALRDR QGNAIGVTAC DIDGDGREEI YFLNTNNAFS GVATYTDKLF KFRNNRWEDI LSDEVNVARG VASLFAGRSV ACVDRKGSGR YSIYIANYAY GNVGPDALIE MDPEASDLR GILALRDVAA EAGVSKYTGG RGVSVGPILS SSASDIFCDN ENGNPFLFHN RGDGTFVDAA ASAGVDDPHQ HGRGVALADF NRDGKVDIVY GNWNGPHRLY LQMSTHGKVR FRDIASPKFS MPSPVRTVIT ADFDNDQELE IFFNNIAYRS SSANRLFRVI RREHGDPLIE ELNPGDALEP EGRGTGGVVT DFDGDGMLDL ILSHGESMAQ PLSVFRGNQG FNNNWLRVVP RTRFGAFARG AKVVLYTKKS GAHLRIIDGG SGYLCEMEPV AHFGLGKDEA SSVEVTWPDG KMVSRNVASG EMNSVLEILY PRDEDTLQDP APLECGQGFS QQENGHCMDT NECIQPFVC PRDKPVCVNT YGSYRCRTNK KCSRGYEPNE DGTACVGTLG QSPGPRPTTP

image



