


Product Name: Recombinant Human PRTFDC1 Protein
Catalog #: PEH4623



Summary

Name	PRDC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRG1
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YGRGVVIMDD WPGYDLNLFT YPQHYYGDLE YVLIPHGIIV DRIERLAKDI MKDIGYSDIM VLCVLKGGYK FCADLVEHLK NISRNSDRFV SMKVDFIRLK SYRNDQSMGE MQIIGDDLS TLAGKNVLIV EDVVTGRTM KALLSNIKEY KPNMIKVASL LVKRTSRSDG FRPDYAGFEI PNLFVVGAL DYNEYFRDLN HICVINEHGK

image



Background

Alternative Names Phosphoribosyltransferase domain-containing protein 1

Note

Product Name: Recombinant Human PRTFDC1 Protein
Catalog #: PEH4623



For Research Use Only , Not for Diagnostic Use