


Product Name: Recombinant Human SDHAF2 Protein
Catalog #: PEH4624



Summary

Name	SDHF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NX18
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SFRRFYRGDS PTDSQKDMIE IPLPPWQERT DESIETKRAR LLYESRKRGML LENCILLSLF AKEHLQHMTE KQLNLYDRLI NEPSNDWDIY YWATEAKPAP EIFENEVMAL LRDFAKNKNK EQRLRAPDLE YLF EKPR

image



Background

Alternative Names	Succinate dehydrogenase assembly factor 2, mitochondrial (SDH assembly factor 2) (SDHAF2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use