



Product Name: Recombinant Human MARCO Protein
Catalog #: PEH4629



Summary

Name	MARCO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UEW3
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVSVRIVGSS NRGRAEVYYS GTWGTICDDE WQNSDAIVFC RMLGYSKGRA LYKVGAGTGQ IWLDNVQCRG TESTLWSCTK NSWGHHDCSH EEDAGVECSV

image



Background

Alternative Names	Macrophage receptor MARCO (Macrophage receptor with collagenous structure) (Scavenger receptor class A member 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use