


**Product Name: Recombinant Human DPP7 Protein**  
**Catalog #: PEH4631**



## Summary

<b>Name</b>	DPP2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & C-His
<b>Accession #</b>	Q9UHL4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DHFNFERFGN KTFPQRFLVS DRFWVRGEGP IFFYTGNEGD VWAFANNSAF VAELAAERGA LLVFAEHRY Y GKSLPFGAQS TQRGHTELLT VEQALADFAE LLRALRRDLG AQDAPAIAFG GSYGGMLSAY LRMKYPHLVA GALAASAPVL AVAGLGDSNQ FFRDVTADFE GQSPKCTQGV REAFRQIKDL FLQGAYDTVR WEFGTCQPLS DEKDLTQLFM FARNAFTVLA MMDYPYPTDF LGPLPANPVK VGCDRLLESEA QRITGLRALA GLVYNASGSE

## image



## Background

<b>Alternative Names</b>	Dipeptidyl peptidase 2 (EC 3.4.14.2) (Dipeptidyl aminopeptidase II) (Dipeptidyl peptidase 7) (Dipeptidyl peptidase II) (DPP II) (Quiescent cell proline dipeptidase)
--------------------------	--

**Product Name: Recombinant Human DPP7 Protein**  
**Catalog #: PEH4631**

---



**Note**

For Research Use Only , Not for Diagnostic Use