


Product Name: Recombinant Human ADGRE2 Protein
Catalog #: PEH4634



Summary

Name	AGRE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UHX3
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGVHSQTLR FFDKVQDLGR DYKPGLANNT IQSILQALDE LLEAPGDLET LPRLQQHCVA SHLLDGLDV LRGLSKNLSN GLLNFSYPAG TELSLEVQKQ VDRSVTLRQN QAVMQLDWNQ AQKSGDPGPS VVGLVSIPGM GKLLAEAPLV LEPEKQMLLH ETHQGLLQDG SPILLSDVIS AFLSNNDTQN

image



Background

Alternative Names	Adhesion G protein-coupled receptor E2 (EGF-like module receptor 2) (EGF-like module-containing mucin-like hormone receptor-like 2) (CD antigen CD312)
--------------------------	--

Note

Product Name: Recombinant Human ADGRE2 Protein
Catalog #: PEH4634



For Research Use Only , Not for Diagnostic Use