


Product Name: Recombinant Human ANGPTL2 Protein
Catalog #: PEH4638



Summary

Name	ANGL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UKU9
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEGSPREFIY LNRYKRAGES QDKCTYTFIV PQRRTGAIC VNSKEPEVLL ENRVHKQELE LLNNEELLKQK RQIETLQQLV EVDGGIVSEV KLLRKESRNM NSRVTQLYMQ LLHEIIRKRD NAEELSQLEN RILNQTADML QLASKYKDLE HKYQHLATLA HNQSEIIAQL EEHCQRVPSA RPVPQPPPAA PPRVYQPPTY NRIINQISTN

image



Background

Alternative Names Angiopoietin-related protein 2 (Angiopoietin-like protein 2)

Note

Product Name: Recombinant Human ANGPTL2 Protein
Catalog #: PEH4638



For Research Use Only , Not for Diagnostic Use