


Product Name: Recombinant Human PLXNB3 Protein
Catalog #: PEH4641



Summary

Name	PLXB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9ULL4
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DAVTGKAKRT LNSRLLRED VEFQPLTMV LVGPGAGGAA GSSEMQRVPA RVLDTDTITQ VKEKVLQVY KGTPFSQRPS VHALDLEWRS GLAGHLTSD EDLTSVTQNH WKRLNTLQHY KVPDGATVGL VPQLHRGSTI SQSLAQRCPL GENIPTLEDG EEGGVCLWHL VKATEEPEGA KVRCSLRER EPARAKAIPE IYLTRLLSMK GTLQKFVDDT FQAILSVMRP IPIAVKYLFD LLELAEKHG

image



Background

Alternative Names Plexin-B3

Note

Product Name: Recombinant Human PLXNB3 Protein
Catalog #: PEH4641



For Research Use Only , Not for Diagnostic Use