


Product Name: Recombinant Human ADAMTS8 Protein
Catalog #: PEH4644



Summary

Name	ATS8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UP79
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRMALYQLDQ QCRQIFGPDF RHCNNTSAQD VCAQLWCHTD GAEPLCHTKN GSLPWADGTP CGPGHLCSEG SCLPEEEVER PKPVADGGWA PWGPWGECSR TCGGGVQFSH RECKDPEPQN GGRYCLGRRA KYQSCHEEC PPDGKSFREQ QCEKYNAYNY

image



Background

Alternative Names	A disintegrin and metalloproteinase with thrombospondin motifs 8 (ADAM-TS 8) (ADAM-TS8) (ADAMTS-8) (EC 3.4.24.-) (METH-2) (METH-8)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use



Product Name: Recombinant Human ADAMTS8 Protein
Catalog #: PEH4644
