


Product Name: Recombinant Human AKT3 Protein
Catalog #: PEH4645



Summary

Name	AKT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y243
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WVQKRGEYIK NWRPRYFLK TDGSFIGYKE KPQDVDLPYP LNNFSVAKCQ LMKTERPKPN TFIIRCLQWT TVIERTFHVD TPEEREEWTE AIQAVADRLQ RQEEERMNCS PTSQIDNIGE EEMDASTTHH

image



Background

Alternative Names RAC-gamma serine/threonine-protein kinase (EC 2.7.11.1) (Protein kinase Akt-3) (Protein kinase B gamma) (PKB gamma) (RAC-PK-gamma) (STK-2)

Note

For Research Use Only , Not for Diagnostic Use