


Product Name: Recombinant Human FIS1 Protein
Catalog #: PEH4646



Summary

Name	FIS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y3D6
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEDLLKFEKK FQSEKAAGSV SKSTQFEYAW CLVRSKYND D IRKGIVLLEE LLPKGSKEEQ RDYVFYLA VG NYRLKEYEKA LKYVRG LLQT EPQNNQAKEL ERLIDKAMKK

image



Background

Alternative Names Mitochondrial fission 1 protein (FIS1 homolog) (hFis1) (Tetratricopeptide repeat protein 11) (TPR repeat protein 11)

Note

For Research Use Only , Not for Diagnostic Use