


Product Name: Recombinant Human ADGRG1 Protein
Catalog #: PEH4647



Summary

Name	AGRG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y653
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLESKLTSVR FMGDMVSFEE DRINATVWKL QPTAGLQDLH IHSRQEEEQS EIMEYSVLLP RTLFQRTKGR SGAEKRLLL VDFSSQALFQ DKNSSQVLGE KVLGIVVQNT

image



Background

Alternative Names	Adhesion G-protein coupled receptor G1 (G-protein coupled receptor 56) (Protein TM7XN1) [Cleaved into: ADGRG1 N-terminal fragment (ADGRG1 NT) (GPR56 N-terminal fragment) (GPR56 NT) (GPR56(N)) (GPR56 extracellular subunit) (GPR56 subunit alpha); ADGRG1 C-terminal fragment (ADGRG1 CT) (GPR56 C-terminal fragment) (GPR56 CT) (GPR56(C)) (GPR56 seven-transmembrane subunit) (GPR56 7TM) (GPR56 subunit beta)]
--------------------------	---

Product Name: Recombinant Human ADGRG1 Protein
Catalog #: PEH4647



Note

For Research Use Only , Not for Diagnostic Use