

Product Name: Recombinant Human SAR1B Protein
Catalog #: PEH4649



Summary

Name	SAR1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6B6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSFIFDWIYS GFSSVLQFLG LYKKTGKLVF LGLDNAGKTT LLHMLKDDRL GQHVPTLHPT SEELTIAGMT FTFDLGGHV QARRVWKNYL PAINGIVFLV DCADHERLLE SKEELDSLMT DETIANVPIL ILGNKIDRPE AISEERLREM FGLYGQTTGK GSISLKELNA RPLEVFMCSV LKRQGYGEGF RWMAQYID

image



Background

Alternative Names GTP-binding protein SAR1b (GTP-binding protein B) (GTBPB)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SAR1B Protein
Catalog #: PEH4649

