


Product Name: Recombinant Human BRAF Protein
Catalog #: PEH4655



Summary

Name	BRAF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15056
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQEEASLAET ALTSGSSPSA PASDSIGPQI LTSPSPSKI PIPQPFRPAD EDHRNQFGQR DRSSAPNVH INTIEPVNID DLIRDQGFRG DGGSTTGLSA TPPASLPGSL TNVKALQKSP GPQREKSSS SSEDNRNMKT

image



Background

Alternative Names Serine/threonine-protein kinase B-raf (EC 2.7.11.1) (Proto-oncogene B-Raf) (p94) (v-Raf murine sarcoma viral oncogene homolog B1)

Note

For Research Use Only , Not for Diagnostic Use