


Product Name: Recombinant Human DNMT1 Protein
Catalog #: PEH4656



Summary

Name	DNMT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P26358
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPDREARAGV QADEDEDGDE KDEKKHRSQP KDLAAKRPE EKEPEKVN PQ ISDEKDEDEK EKRRTTPK EPTEKKMARA KTMNSKTHP PKCIQCGQYL DDPDLKYGQH PPDVDEPQM LTNEKLSIFD ANESGFESYE ALPQHKLTCF SVYCKHGHLC PIDTGLIEKN

image



Background

Alternative Names	DNA (cytosine-5)-methyltransferase 1 (Dnmt1) (EC 2.1.1.37) (CXXC-type zinc finger protein 9) (DNA methyltransferase Hsa1) (DNA MTase Hsa1) (M.Hsa1) (MCMT)
--------------------------	--

Note

Product Name: Recombinant Human DNMT1 Protein
Catalog #: PEH4656



For Research Use Only , Not for Diagnostic Use