



Product Name: Recombinant Human KRT7 Protein
Catalog #: PEH4661



Summary

Name	K2C7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08729
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFASFIDKVR FLEQQNKLLE TKWTLLQEQK SAKSSRLPDI FEAQIAGLRG QLEALQVDGG RLEAELRSMQ DVVEDFKNKY EDEINHRTAA ENEFVVLKGD VDAAYMSKVE LEAKVDALND EINFLRTLNE TELTELQSQI SDTSVVLMSD NSRSLDLGI IAEVKAQYEE MAKCSRAEAE AWYQTKFETL QAQAGKHGDD LRNTRNEISE MNRAIQRLQA EIDNIKNQRA KLEAAIAEAE ERGELALKDA RAKQEELEAA LQRGKQDMAR QLREYQELMS VKLALDIEIA TYRKLLEGEE SRLAGDGVGA VNISVMNSTG

image



Background

Alternative Names Keratin, type II cytoskeletal 7 (Cytokeratin-7) (CK-7) (Keratin-7) (K7) (Sarcolectin) (Type-II keratin Kb7)

Product Name: Recombinant Human KRT7 Protein
Catalog #: PEH4661



Note

For Research Use Only , Not for Diagnostic Use