


Product Name: Recombinant Human ARM CX6 Protein
Catalog #: PEH4662



Summary

Name	ARMX6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L4S7
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEEEWDDDDQE LDEEPPDIWF DFETMARPWT EDGDWTEPGA PGGTEDRPSG GGKANRAHPI KQRPFYEHK NTWSAQNCKN GSCVLDLSKC LFIQGKLLFA EPKDAGFPFS QDINSHLASL SMARNTSPTP DPTVREALCA PDNLNASIES

image



Background

Alternative Names Protein ARM CX6

Note

For Research Use Only , Not for Diagnostic Use