


Product Name: Recombinant Human EPOP Protein
Catalog #: PEH4675



Summary

Name	EPOP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NHQ4
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AACPRRGEEE EGGGGFPHFG VRSCAPPGRC PAPHPREST TSFASAPPRP APGLEPQRGP AASPPQEPSS RPPSPAGLS TEPAGPGTAP RPFLPGQPAE VDGNPPPAAP EPAASPSTA SPAPAAPGDL RQEHFDRLIR RSKLWCYAKG FALDTPSLRR GPERPPAKGP ARGAAKKRRL PAPPPRTAQP

image



Background

Alternative Names	Elongin BC and Polycomb repressive complex 2-associated protein (Proline-rich protein 28)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EPOP Protein
Catalog #: PEH4675

