


Product Name: Recombinant Human SPATA31C2 Protein
Catalog #: PEH4691



Summary

Name	S31C2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	B4DYI2
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LMEEAVSEFE PGKATKSETQ PQVSATVVLL PDGQASVVPH ASENLASQVP QGHLLQSMPTG NMQASQELCD LMSARRSNMG HKEPRNPNCQ GSCKSQSPMF PPTHKRENSR KPNLEKHEEM FQGLRTPQLT PGRKTEDTRQ NEGVQLLPSK KQPPSISHFG

image



Background

Alternative Names Putative spermatogenesis-associated protein 31C2 (Protein FAM75C2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SPATA31C2 Protein
Catalog #: PEH4691

