


Product Name: Recombinant Human USP17L15 Protein
Catalog #: PEH4693



Summary

Name	U17LF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	C9J2P7
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRATQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTKEPF NVRKVEGTLP PDVLVIHQSK YKCGMKNHHP EQQSSLLNLS STTPTHQESM NTGTLASLRG RARRSKGKNK HSKRALLVCQ

image



Background

Alternative Names Ubiquitin carboxyl-terminal hydrolase 17-like protein 15 (EC 3.4.19.12)

Note

For Research Use Only , Not for Diagnostic Use