


Summary

Name	MEIOS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	C9JSJ3
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVEGEDKVNP SEPHGLRMEE KWLLKGKLRN QRNQNKLLSP NKKQRKNHTS KLQELALLP IALKTGTKKL TKKEILVHVL QYIQYLQRNI DAAKALFKCH ITTGEGGLAG LGQKPAWGPA RRRRHSTPSS SPSSQKSLQ GACQKPRKKK LTQASESQTR TPKPRSLAL NKPEKLVAPS

image



Background

Alternative Names	Meiosis initiator protein (Basic helix-loop-helix and HMG box domain-containing protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MEIOSIN Protein
Catalog #: PEH4696

