


Product Name: Recombinant Human USP17L18 Protein
Catalog #: PEH4699



Summary

Name	U17LI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	D6R9N7
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRALGAEDTD RRAKQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTKPEF NVRKVEGTLP PDVLVIHQSK YKCGMKNHHP EQQSLLNLS STTPHQESM NTGTLASLRG RARRSKGKNK HSKRALLVCQ

image



Background

Alternative Names Ubiquitin carboxyl-terminal hydrolase 17-like protein 18 (EC 3.4.19.12)

Note

For Research Use Only , Not for Diagnostic Use