

Product Name: Recombinant Human TAF4 Protein
Catalog #: PEH4705



Summary

Name	TAF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00268
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEPGGGSFRD DDDINDVASM AGVNLSEESA RILATNSELV GTLTRSCKDE TFLQAPLQR RILEIGKKHG ITLHPDVVS YVSHATQQRL QNLVEKISSET AQQKNFSYKD DDRYEQASDV RAQLKFFEQL DQIEKQRKDE QEREILMRAA KRSRQEDPE QLRLKQKAKE MQQQELAQMR

image



Background

Alternative Names	Transcription initiation factor TFIID subunit 4 (RNA polymerase II TBP-associated factor subunit C) (TBP-associated factor 4) (Transcription initiation factor TFIID 130 kDa subunit) (TAF(II)130) (TAFII-130) (TAFII130) (Transcription initiation factor TFIID 135 kDa subunit) (TAF(II)135) (TAFII-135) (TAFII135)
--------------------------	---

Product Name: Recombinant Human TAF4 Protein
Catalog #: PEH4705



Note

For Research Use Only , Not for Diagnostic Use