


Summary

Name	DLGP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14490
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESLDSMKALT AAIEAANAQI HGPASQHMGN NTATVTTTTT IATVTTEDRK KDHFKKNRCL SIGIQVDDAE EPDKTGENKA PSKFQSVGQVQ VEEKCFRRF TRNSVTTAV QADLDFHDNL ENSLESIEDN SCPGPMARQF SRDASTSTVS IQGSGNHYHA CAADDDFDTD

image



Background

Alternative Names	Disks large-associated protein 1 (DAP-1) (Guanylate kinase-associated protein) (hGKAP) (PSD-95/SAP90-binding protein 1) (SAP90/PSD-95-associated protein 1) (SAPAP1)
--------------------------	--

Note

Product Name: Recombinant Human DLGAP1 Protein
Catalog #: PEH4709



For Research Use Only , Not for Diagnostic Use