


Product Name: Recombinant Human NRG2 Protein
Catalog #: PEH4710



Summary

Name	NRG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14511
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKKEVGKILC TDCATRPKLIK KMKSTGQVG EKQSLKCEAA AGNPQPSYRW FKDGGKELNRS RDIRIKYGNG RKNSRLQFNK VKVEDAGEYV CEAENILGKD TVRGRLYVNS

image



Background

Alternative Names Pro-neuregulin-2, membrane-bound isoform (Pro-NRG2) [Cleaved into: Neuregulin-2 (NRG-2) (Divergent of neuregulin-1) (DON-1) (Neural- and thymus-derived activator for ERBB kinases) (NTAK)]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NRG2 Protein
Catalog #: PEH4710

