


Product Name: Recombinant Human UTY Protein
Catalog #: PEH4715



Summary

Name	UTY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14607
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLTLPSNSVP QGDADSHLSC HTATSGGQQG IMFTKESKPS KNRSVPETS RHTGDTSNGC ADVKGLSNHV HQLIADAVSS PNHGDSPNLL IADNPQLSAL LIGKANGNVG TGTCDKVNNI HPAVHTKTDH SVASSPSSAI STATPSPKST EQRSINSVTS LNSPHSGLHT VNGEGLGKSQ SSTKVDLPLA SHRSTSQILP SMSVSICPSS TEVLKACRNP

image



Background

Alternative Names	Histone demethylase UTY (EC 1.14.11.68) (Ubiquitously-transcribed TPR protein on the Y chromosome) (Ubiquitously-transcribed Y chromosome tetratricopeptide repeat protein) ([histone H3]-trimethyl-L-lysine(27) demethylase UTY)
--------------------------	---

Product Name: Recombinant Human UTY Protein
Catalog #: PEH4715



Note

For Research Use Only , Not for Diagnostic Use