


Product Name: Recombinant Human BTAF1 Protein
Catalog #: PEH4716



Summary

Name	BTAF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14981
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKSGEVDPK E RIARQRKLLQ KKLGLNMGEA IGMSTEELFN DEDLDYTPTS ASFVNKQPTL QAAELIDSEF RAGMSNRQKN KAKRMAKLFA KQRSRDVET NEKSNDSTDG EPEEKRRKIA NVVINQSAND SKVLIDNIPD SSSLIEETNE WPLESFCEEL

image



Background

Alternative Names	TATA-binding protein-associated factor 172 (EC 3.6.4.-) (ATP-dependent helicase BTAF1) (B-TFIID transcription factor-associated 170 kDa subunit) (TAF(II)170) (TBP-associated factor 172) (TAF-172)
--------------------------	---

Note

Product Name: Recombinant Human BTAF1 Protein
Catalog #: PEH4716



For Research Use Only , Not for Diagnostic Use