



Product Name: Recombinant Human U2SURP Protein
Catalog #: PEH4717



Summary

Name	SR140_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15042
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDEEDTQSSK SEEHLYSNP IKEEMTESKF SKYSEMSEEK RAKLREIELK VMKFQDELES GKRPKKPGQS FQEQVEHYRD KLLQREKEKE LERERERDKK DKEKLESRSK DKKEKDECTP TRKERKRRHS TSPSPSRSSS GRRVKSPSPK SERSERSERS

image



Background

Alternative Names U2 snRNP-associated SURP motif-containing protein (140 kDa Ser/Arg-rich domain protein) (U2-associated protein SR140)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human U2SURP Protein
Catalog #: PEH4717

