


Product Name: Recombinant Human KDM6A Protein
Catalog #: PEH4720



Summary

Name	KDM6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15550
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAGPNGERPL SSTGPSQHLQ AAGSGIQNQNH GHPTLPSNSV TQGAALNHLS SHTATSGGQQ GITLTKEKSP SGNILTVPET SRHTGETPNS TASVEGLPNH VHQMATAVC SPSHGDSKSP GLLSSDNPQL SALLMGKANN NVGTGTCDKV NNIHPAVHTK TDNSVASSPS SAISTATPSP KSTEQTTTNS VTSLNSPHSG

image



Background

Alternative Names	Lysine-specific demethylase 6A (EC 1.14.11.68) (Histone demethylase UTX) (Ubiquitously-transcribed TPR protein on the X chromosome) (Ubiquitously-transcribed X chromosome tetratricopeptide repeat protein) ([histone H3]-trimethyl-L-lysine(27) demethylase 6A)
--------------------------	---

Product Name: Recombinant Human KDM6A Protein
Catalog #: PEH4720



Note

For Research Use Only , Not for Diagnostic Use