


Product Name: Recombinant Human MGAM Protein
Catalog #: PEH4724



Summary

Name	MGA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43451
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGELFWDNGE TKDTVANKVY LLCEFSVTQN RLEVNISQST YKDPNNLAFN EIKILGTEEP SNVTVKHNGV PSQTSPTVTY DSNLKVAIIT DIDLLLGEAY TVEWSIKIRD EEKIDCYPDE NGASAENCTA RGCWEASNS SGVPFCYFVN DLYSVSDVQY

image



Background

Alternative Names Maltase-glucoamylase, intestinal [Includes: Maltase (EC 3.2.1.20) (Alpha-glucosidase); Glucoamylase (EC 3.2.1.3) (Glucan 1,4-alpha-glucosidase)]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MGAM Protein
Catalog #: PEH4724

