


Product Name: Recombinant Human SMARCA5 Protein
Catalog #: PEH4730



Summary

Name	SMCA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60264
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VQDFQFFPPR LFELEKEIL FYRKTIGYKV PRNPELPNAA QAQKEEQLKI DEAESLNDEE LEEKEKLLTQ GFTNWNKRDF NQFIKANEKW GRDDIENIAR EVEGKTPEEV IEYSAVFWER CNEHQDIEKI

image



Background

Alternative Names	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 5 (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin A5) (EC 3.6.4.-) (Sucrose nonfermenting protein 2 homolog) (hSNF2H)
--------------------------	---

Note

Product Name: Recombinant Human SMARCA5 Protein
Catalog #: PEH4730



For Research Use Only , Not for Diagnostic Use