


Product Name: Recombinant Human ST18 Protein
Catalog #: PEH4731



Summary

Name	ST18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60284
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NNVFVTHRSL SGCPLNAQVI KKGKVSEELM TIKLKATGGI ESDEEIRHLD EEIKELNESN LKIEADMMKL QTQITSMESN LKTIEEENKL IEQNNESSLK ELAGLSQALI SSLADIQLPQ MGPISEQNFE AYVNTLTDMY SNLERDYSPE

image



Background

Alternative Names Suppression of tumorigenicity 18 protein (Zinc finger protein 387)

Note

For Research Use Only , Not for Diagnostic Use