


Product Name: Recombinant Human OPA1 Protein
Catalog #: PEH4734



Summary

Name	OPA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60313
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEHDDIFDKL KEAVKEESIK RHKWNDFAED SLRVIQHNAL EDRSISDKQQ WDAAIYFMEE ALQARLKDTE NAIENMVGPD WKKRWLYWKN RTQEQCVHNE TKNELEKMLK

image



Background

Alternative Names	Dynamin-like 120 kDa protein, mitochondrial (EC 3.6.5.5) (Optic atrophy protein 1) [Cleaved into: Dynamin-like 120 kDa protein, form S1]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use