


Product Name: Recombinant Human TRIM13 Protein
Catalog #: PEH4738



Summary

Name	TRI13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60858
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFESLFQSFE TWRRGDALSR LDTLETSKRK SLQLLTKDSD KVKEFFEKLQ HTLDQKKNEI LSDFETMKLA VMQAYDPEIN KLNTILQEQR MAFNIAEAFK DVSEPIVFLQ QMQEFREKIK VIKETPLPPS NLPASPLMKN

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM13 (EC 2.3.2.27) (B-cell chronic lymphocytic leukemia tumor suppressor Leu5) (Leukemia-associated protein 5) (Putative tumor suppressor RFP2) (RING finger protein 77) (RING-type E3 ubiquitin transferase TRIM13) (Ret finger protein 2) (Tripartite motif-containing protein 13)
--------------------------	--

Product Name: Recombinant Human TRIM13 Protein
Catalog #: PEH4738



Note

For Research Use Only , Not for Diagnostic Use