

Summary

Name	LIRB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75022
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFQRPAGAAE TEPKDRGLLR RSPPAADVQE ENLYAAVKDT QSEDRVELDS QSPHDEDPQA VTYAPVKHSS PRREMASPPS SLSGFELDTK DRQVEEDRQM DTEAAASEAS QDVTYAQLHS LTLRRKATEP PPSQEGEPPA EPSIYATLAI

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3) (Leukocyte immunoglobulin-like receptor 3) (CD85 antigen-like family member A) (Immunoglobulin-like transcript 5) (ILT-5) (Monocyte inhibitory receptor HL9) (CD antigen CD85a)
--------------------------	---

Note

Product Name: Recombinant Human LILRB3 Protein
Catalog #: PEH4739



For Research Use Only , Not for Diagnostic Use