


Product Name: Recombinant Human ZNF217 Protein
Catalog #: PEH4743



Summary

| | |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | ZN217_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | O75362 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 25 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | EGIHLDKNDD GGIKHLTSS RECSYCGKFF RSNYYLNIHL RTHTGEKPYK CEFCEYAAAQ KTSRLRYHLER HHKEKQTDVA AEVKNDGKNQ DTEDALLTAD SAQTKNLKRF FDGAKDVTGS PPAKQLKEMP SVFQNVLGSA VLSPAHKDTQ DFHKNAADDS ADKVNKNPTP AYLDLLKKRS AVETQANNLI CRTKADVTTP PDGSTTHNLE |

image



Background

Alternative Names Zinc finger protein 217

Note

Product Name: Recombinant Human ZNF217 Protein
Catalog #: PEH4743



For Research Use Only , Not for Diagnostic Use