


Product Name: Recombinant Human SPAG7 Protein
Catalog #: PEH4744



Summary

Name	SPAG7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75391
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SMEKPPSLGD QETRRKAREQ AARLKKLQEQ EKQQKVEFRK RMEKEVSDFI QDSGQIKKKF QPMNKIERSI LHDVVEVAGL TSFSFGEDDD CRYVMIFKKE FAPSDEELDS YRRGEEWDPQ KAEKRKLKE LAQRQEEEA QGPPVVVSPA

image



Background

Alternative Names Sperm-associated antigen 7

Note

For Research Use Only , Not for Diagnostic Use