


Product Name: Recombinant Human POLQ Protein
Catalog #: PEH4745



Summary

Name	DPOLQ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75417
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLIKKSNVNE NQDTHQQLTC SNDESIIFSE MDSVQMVEAL DNVDIFPVQE KNHTVVSPRA LELSDPVLDE HHQGDQDGGD QDERAEKSKL TGTRQNHFSI WSGASFDLSP GLQRILDKVS SPLENEKLKS MTINFSSLNR KNTLNREEQE VISNLETKQV QGISFSSNNE

image



Background

Alternative Names DNA polymerase theta (EC 2.7.7.7) (DNA polymerase eta)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POLQ Protein
Catalog #: PEH4745

