


Product Name: Recombinant Human ZNF324 Protein
Catalog #: PEH4747



Summary

Name	Z324A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75467
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPERQKPCAQ EVPGRTEFGSA QDLEAAGGRG HHRMGAVWQE PHRLGGQEP STWDELGEAL HAGEKSFECR ACSKVFKSS DLLKHLRHT GERPYECAQC GKAFSQTSHL

image



Background

Alternative Names Zinc finger protein 324A (Zinc finger protein ZF5128)

Note

For Research Use Only , Not for Diagnostic Use