


Product Name: Recombinant Human ZNF821 Protein
Catalog #: PEH4749



Summary

Name	ZN821_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75541
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTDFPGDLGS QRQAIQQLRD QDSSSSDSEG DEEETTQDEV SSHTSEEDGG VVKVEKELEN TEQPVGGNEV VEHEVTGNLN SDPLLELCQC PLCQLDCGSR EQLIAHVYQH TAAVVSASY MCPVCGRALS SPGSLGRHLL IHSEDQRSNC AVCGARFTSH ATFNSEKLPE VLNMESLPTV HNEGPPSAEG KDIAFSPPVY PAGILLVCNN CAAYRKLLEA QTPSVRKWAL RRQNEPLEVR LQRLERERTA KKSRRDNETP

image



Background

Alternative Names Zinc finger protein 821

Product Name: Recombinant Human ZNF821 Protein
Catalog #: PEH4749



Note

For Research Use Only , Not for Diagnostic Use