


Product Name: Recombinant Human WDHD1 Protein
Catalog #: PEH4752



Summary

Name	WDHD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75717
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSVIFHNHLD YLAKNGYEYE ESTKNQATKE QQELLMKMLA LSCKLEREFR CVELADLMTQ NAVNLAIKYA SRSRKLILAQ KLSELAVEKA AELTATQVEE EEEEEDFRKK LNAGYSNTAT EWSQPRFRNQ VEEDAEDSGE ADDEEKPEIH KPGQNSFSKS TNSSDVSAKS GAVTFSSQGR VNPFKVSASS KEPAMSMNSA

image



Background

Alternative Names	WD repeat and HMG-box DNA-binding protein 1 (Acidic nucleoplasmic DNA-binding protein 1) (And-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human WDHD1 Protein
Catalog #: PEH4752

