


Product Name: Recombinant Human KRT37 Protein
Catalog #: PEH4758



Summary

Name	KRT37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O76014
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNDRLANYLE KVRQLEQENA ELETLLERS KCHESTVCPD YQSYFRTIEE LQQKILCSKA ENARLIVQID NAKLAADDFR IKLESERSLH QLVEADKCGT QKLLDDATLA KADLEAQQES LKEEQLSLKS NHEQEVKILR SQLGEKFRIE LDIEPTIDLN RVLGEMRAQY EAMVETNHQD VEQWFQAQSE GISLQAMSCS EELQCCQSEI LELRCTVNAL EVERQAQHTL KDCLQNSLCE AEDRYGTELA QMQSLISNLE EQLSEIRADL ERQNQEYQVL LDVKARLENE IATYRNLLES

image



Background

Alternative Names	Keratin, type I cuticular Ha7 (Hair keratin, type I Ha7) (Keratin-37) (K37)
--------------------------	---

Product Name: Recombinant Human KRT37 Protein
Catalog #: PEH4758



Note

For Research Use Only , Not for Diagnostic Use