


Product Name: Recombinant Human KRT38 Protein
Catalog #: PEH4759



Summary

Name	KRT38_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O76015
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGHEKETMQF LNDRLANYLE KVRQLEQENA ELEATLLERS KCHESTVCPD YQSYFHTIEE LQQKILCSKA ENARLIVQID NAKLAADDFR IKLESERSLR QLVEADKCGT QKLLDDATLA KADLEAQQES LKEEQLSLKS NHEQEVKILR SQLGEKLRIE LDIEPTIDLN RVLGEMRAQY EAMLETNRQD VEQWFQAQSE GISLQDMSCS EELQCCQSEI LELRCTVNAL EVERQAQHTL KDCLQNSLCE AEDRFGTELA QMQSLISNVE EQLSEIRADL ERQNQEYQVL LDVKTRLENE IATYRNLLS EDCKLPCNPC

image



Background

Alternative Names	Keratin, type I cuticular Ha8 (Hair keratin, type I Ha8) (Keratin-38) (K38)
--------------------------	---

Product Name: Recombinant Human KRT38 Protein
Catalog #: PEH4759



Note

For Research Use Only , Not for Diagnostic Use