


**Product Name: Recombinant Human UBR5 Protein**  
**Catalog #: PEH4767**



## Summary

<b>Name</b>	UBR5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95071
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	19 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EHENDDDTNQ SATLNDKDDD SLPAETGQNH PFFRRSDSMT FLGCIPPNPF EVPLAEAIPL ADQPHLLQPN ARKEDLFGRP SQGLYSSAS SGKCLMEVTV DRNCLEVLPT KMSYAANLKN VMNMQRQKK EGEEQVLP ETESSKPGPS

## image



## Background

**Alternative Names** E3 ubiquitin-protein ligase UBR5 (EC 2.3.2.26) (E3 ubiquitin-protein ligase, HECT domain-containing 1) (HECT-type E3 ubiquitin transferase UBR5) (Hyperplastic discs protein homolog) (hHYD) (Progesterin-induced protein)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human UBR5 Protein**  
**Catalog #: PEH4767**

---

