


Product Name: Recombinant Human SCAF4 Protein
Catalog #: PEH4768



Summary

Name	SCAF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95104
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRHRDLEERN RRSSGHRDRE RDSRDRESRR EKEEARGKEK PEVTDRAAGN KTVEPPISQV GNVDTASELE KGVSEAAVLK

image



Background

Alternative Names	SR-related and CTD-associated factor 4 (CTD-binding SR-like protein RA4) (Splicing factor, arginine/serine-rich 15)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use