


Product Name: Recombinant Human MBD3 Protein
Catalog #: PEH4773



Summary

Name	MBD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95983
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDQPRQLFWE KKL SGLNAFD IAEELVKTMD LPKGLQGVGP GCTDETL LSA IASALHTSTM PITGQLSAAV EKNPGVWLNT TQPLCKAFMV TDEDIRKQEE LVQQVRKRLE EALMADMLAH VEELARDGEA PLDKACAEDD DEEDEEEEE

image



Background

Alternative Names Methyl-CpG-binding domain protein 3 (Methyl-CpG-binding protein MBD3)

Note

For Research Use Only , Not for Diagnostic Use